

F61: Nuclear Magnetic Resonance

T. Gierlich und A. Impertro

May 26th, 2017

Outline

Introduction and Theoretical Concepts

Part I: Relaxation Times

Introduction

- ▶ **Nuclear:** Interaction of the nuclear spin...
- ▶ **Magnetic:** ...with magnetic fields
- ▶ **Resonance:** → resonant interaction
- ▶ **Spectroscopy:** Resolve the different signal components

Applications

- ▶ Non-destructive detection of substances in a sample
- ▶ Determination of molecular structures and dynamics
- ▶ Multidimensional imaging

Theoretical Concepts - Working Principle




Figure: Magnetization

- ▶ Nuclei with spin J have a magnetic moment μ
- ▶ Minimal energy \rightarrow Dipole aligned parallel to B-field
- ▶ Ensemble of many nuclei: Measurable magnetization \vec{M}
- ▶ Ground state $\rightarrow M_{\perp} = 0$




Figure: Larmor-Precession of M_{\perp}

- ▶ Excited states have a component $M_{\perp} \neq 0$
- ▶ M_{\perp} precesses around the field lines with the Larmor frequency

$$\omega_L = \gamma B_0 \quad (1)$$

- ▶ ω_L can be measured!

Theoretical Concepts - Working Principle

How can we create an excited state ?

- ▶ An oscillating B-Field \vec{B}_1 rotates the magnetization \vec{M} by an angle

$$\alpha = \gamma B_1 \Delta t \quad (2)$$

- ▶ By choosing Δt , we can create:
 - ▶ A perpendicular magnetization (90°-Pulse)
 - ▶ An anti-parallel magnetization (180°-Pulse)




Figure: Rotation of M due to a HF-Pulse

Theory of Relaxation

Excited states decay into the Ground State on a characteristic timescale.
The decay is of exponential nature and described in the *Bloch equations*:

$$\frac{dM_{\perp}(t)}{dt} = -\frac{M_{\perp}(t)}{T_2} \quad (3)$$

$$\frac{dM_{\parallel}(t)}{dt} = -\frac{M_{\parallel}(t) - M_0}{T_1} \quad (4)$$

- ▶ T_2 : Spin-Lattice Relaxation
- ▶ T_1 : Spin-Spin Relaxation

T_2 -Measurement: Spin Echo

Spin-Echo principle




Figure: Principle of the spin-echo method

Pulse sequence




Figure: Spin-Echo measurement with $\tau = 10\text{ms}$

- **Disadvantage:** Dephasing for long measurement times!

T_1 -Measurement: Spin Echo




Figure: T2-Messung Probe 1 mit Fit.




Figure: T2-Messung Probe 3 mit Fit.

T_2 -Measurement: Carr-Purcell Sequence

Improve dephasing problem of spin-echo method:

- ▶ Inject a 180° -Pulse on odd multiples of a time τ .
- ▶ The system is phase coherent on even multiples of a time τ .




Figure: T2-Messung über Carr-Purcell Probe 1 mit Fit.




Figure: T2-Messung über Carr-Purcell Probe 3 mit Fit.

T_1 -Measurement: Spin Echo

Spin-Echo, but start with a 180° -Pulse (Anti-parallel Magnetization)




Figure: T1-Messung Probe 1 mit Fit.




Figure: T1-Messung Probe 3 mit Fit.

Relaxation Times: Evaluation

Table: Relaxation times: Measured values

Zeit Methode	T_1 [ms] Spin-Echo	T_2 [ms] Spin-Echo	T_2 [ms] Carr-Purcell
Probe 1 (Gd 1:500)	$(125,5 \pm 0,6)$	$(119,5 \pm 0,5)$	$(140,1 \pm 0,4)$
Probe 3 (Gd 1:600)	$(150,5 \pm 1,2)$	$(139,3 \pm 0,8)$	$(166,9 \pm 0,4)$