

F61: NMR–Spectroscopy

T. Gierlich and A. Impertro

May 26th, 2017

Outline

Introduction and Theoretical Concepts

Part I: Relaxation Times

Introduction

NMR-Spectroscopy

Nuclear Magnetic Resonance

The diagram illustrates the components of NMR-Spectroscopy. The main title 'NMR-Spectroscopy' is at the top. Three arrows point from the letters 'N', 'M', and 'R' respectively to the words 'Nuclear', 'Magnetic', and 'Resonance' below. The words are in a bold, sans-serif font.

Introduction

NMR-Spectroscopy

Nuclear

Magnetic

Resonance

Detection of substances




Figure: Carl Nave, Hyperphysics, hyperphysics.phy-astr.gsu.edu

Multidimensional imaging




Figure: Sierra Vista Diagnostics, svdrads.com

Theoretical Concepts - Working Principle




Figure: Magnetization

- ▶ Nuclei with spin I have a magnetic moment μ
- ▶ Ensemble of many nuclei: Measurable magnetization \vec{M}
- ▶ Minimal energy \rightarrow Dipole aligned parallel to B-field
- ▶ Ground state $\rightarrow M_{\perp} = 0$

Theoretical Concepts - Working Principle




Figure: Magnetization

- ▶ Nuclei with spin I have a magnetic moment μ
- ▶ Ensemble of many nuclei: Measurable magnetization \vec{M}
- ▶ Minimal energy \rightarrow Dipole aligned parallel to B-field
- ▶ Ground state $\rightarrow M_{\perp} = 0$




Figure: Larmor-Precession of M_{\perp}

- ▶ Excited states have a component $M_{\perp} \neq 0$
- ▶ M_{\perp} precesses around the field lines with the Larmor frequency

$$\omega_L = \gamma B_0 \quad (1)$$

- ▶ ω_L can be measured!

Theoretical Concepts - Working Principle

How can we create an excited state ?

- ▶ An oscillating B-Field \vec{B}_1 rotates the magnetization \vec{M} by an angle

$$\alpha = \gamma B_1 \Delta t \quad (2)$$




Figure: Rotation of M due to an HF-Pulse

Theoretical Concepts - Working Principle

How can we create an excited state ?

- ▶ An oscillating B-Field \vec{B}_1 rotates the magnetization \vec{M} by an angle

$$\alpha = \gamma B_1 \Delta t \quad (2)$$




Figure: Rotation of M due to an HF-Pulse

- ▶ By choosing Δt , we can create:
 - ▶ A perpendicular magnetization (90°-Pulse)
 - ▶ An anti-parallel magnetization (180°-Pulse)

Setup and Measurement Principle




Figure: McGraw Hill Higher Education, mhhe.com

Theory of Relaxation

Excited states decay into the Ground State on a characteristic timescale.
The decay is of exponential nature and described in the *Bloch equations*:

$$\frac{dM_{\perp}(t)}{dt} = -\frac{M_{\perp}(t)}{T_2} \quad (3)$$

$$\frac{dM_{\parallel}(t)}{dt} = -\frac{M_{\parallel}(t) - M_0}{T_1} \quad (4)$$

- ▶ T_2 : Spin-Lattice Relaxation
- ▶ T_1 : Spin-Spin Relaxation

T_2 -Measurement: Spin Echo

Spin-Echo principle




Figure: Principle of the spin-echo method

T_2 -Measurement: Spin Echo

Spin-Echo principle




Figure: Principle of the spin-echo method

Pulse sequence




Figure: Spin-Echo measurement with $\tau = 10\text{ms}$

T_2 -Measurement: Spin Echo

Spin-Echo principle




Figure: Principle of the spin-echo method

Pulse sequence




Figure: Spin-Echo measurement with $\tau = 10\text{ms}$

- **Disadvantage:** Dephasing for long measurement times!

T_2 -Measurement: Spin Echo




Figure: T2-Measurement Sample 1 with fit.




Figure: T2-Measurement Sample 3 with fit.

T_2 -Measurement: Carr-Purcell Sequence

Improve dephasing problem of spin-echo method:

- ▶ Inject a 180° -Pulse on odd multiples of a time τ .
- ▶ The system is phase coherent on even multiples of a time τ .

T_2 -Measurement: Carr-Purcell Sequence

Improve dephasing problem of spin-echo method:

- ▶ Inject a 180° -Pulse on odd multiples of a time τ .
- ▶ The system is phase coherent on even multiples of a time τ .




Figure: T2-Measurement using Carr-Purcell, Sample 1, with fit.




Figure: T2-Measurement using Carr-Purcell, Sample 3, with fit.

T_1 -Measurement

Start with a 180° -Pulse (Anti-parallel Magnetization) and probe the magnetization after time τ with a 90° -Pulse

T_1 -Measurement

Start with a 180° -Pulse (Anti-parallel Magnetization) and probe the magnetization after time τ with a 90° -Pulse




Figure: T1-Measurement Sample 1 with fit.




Figure: T1-Measurement Sample 3 with fit.

Relaxation Times: Evaluation

Table: Relaxation times – Measured values

Time Method	T_1 [ms] 180°-90°	T_2 [ms] Spin-Echo	T_2 [ms] Carr-Purcell
Sample 1 (Gd 1:500)	(125, 5 ± 0, 6)	(119, 5 ± 0, 5)	(140, 1 ± 0, 4)
Sample 3 (Gd 1:600)	(150, 5 ± 1, 2)	(139, 3 ± 0, 8)	(166, 9 ± 0, 4)